

FLOCATH *quick*™ Intermittent Catheter

Your Solution for Quick and Easy Intermittent Catheterization!

Hydrophilic-Coated Catheter for Self-Catheterization

The Rüscher FLOCATH Quick™ is an "All-in-One" application providing the ultimate convenience for intermittent catheterization. The integrated package contains a hydrophilic coated flexible PVC catheter and a sterile saline solution pouch. Our patented sheath maintains and protects the sterile catheter during insertion.

The FLOCATH Quick™ with its hydrophilic coating glides easily into the urethra after being activated by the 0.9% sterile saline solution. The unique protection sheath facilitates sterile catheterization technique so no additional gloves or messy gels are needed.




FLOCATH Quick™ comes pre-packaged in sterile saline solution and in an easy-to-open foil container that can be taken anywhere.

Easy opening/handling tab
Allows the packaging to be opened more easily, even when manual dexterity is limited. It can also be used to hang up the catheter pack.

Self-adhesive patch
Located on the back of the package, it enables the catheter to be stuck onto a suitable surface, allowing both hands to remain free.

Integrated sterile saline solution pouch
Whether at work, rest or play, sterile saline solution is ready at hand to activate the hydrophilic coating without any hassle or mess.

FLOCATH <i>quick</i> ™		
Product #	Size	Quantity
220400060	6 Fr	30
220400080	8 Fr	30
220400100	10 Fr	30
220400120	12 Fr	30
220400140	14 Fr	30
220400160	16 Fr	30




Simply release the sterile saline from its pouch and the hydrophilic coating will be activated after 30 seconds.

Water binds to the entire hydrophilic coated area, creating a protective, slippery smooth layer of water.

Ultra-smooth catheter eyelets
Provide maximum comfort on insertion and removal. FLOCATH Quick™ provides 90% less friction during insertion/withdrawal compared to a traditional gel-lubricated catheter.

Our unique patented protective catheter sleeve
The catheter can be inserted safely and hygienically without directly touching it. After insertion, the protective sleeve can be pushed back and used as an extension to the catheter.



Latex-Free

To place an order contact
Customer Service at **1-800-553-5214**

Rüscher Inc • 2450 Meadowbrook Parkway
Duluth GA • 30096 • www.rusch.com
P: 770 623 0816 • F: 770 623 1829

RÜSCH
A TELEFLEX COMPANY