




# 2XL BRIEF PRODUCT USE STANDARD GUIDELINE

The purpose of this guide is to provide a standard for facilities for the use of the First Quality 2XL brief: **PV-017**.

- The 2XL brief is a specialty product designed for a specific population within the facility.
- It is not to be used on the general resident population.
- Prior to placing the 2XL brief on a resident the below criteria should be met.

## 2XL BRIEF CRITERIA/INDICATIONS:

- The waist measurement of the residents wearing this brief are between **62"-73"**
- The typical weight of the residents wearing this brief is greater than **260-300 lbs.**
- The 2X size is not comparable to clothing you would purchase that is labeled 2X, this 2X size is a double extra large adult brief.
- In the event the resident cannot be weighed, utilize the First Quality tape measure specific for 2XL briefs.
- The 2XL brief should not be stored with the other briefs; it should be placed directly in the resident's room.
- Criteria for ordering and using the 2XL Brief:
  - The First Quality XL brief has been tried and does not fit the resident
  - The resident is unable to utilize regular underwear with a pad or liner.
  - The Director of Nursing at the facility has approved the use of the 2XL Brief and communicated approval with Central Supply.



Height & Weight Chart

## FIRST QUALITY SIZING TOOLS

First Quality's Height & Weight Sizing Charts cover sizes Small to XLarge and up to 260 lbs.

If your resident exceeds the parameters of the Height and Weight Chart, please use First Quality's Sizing Tape Measure. Residents who exceed 260 lbs should be measured hip to hip for accurate sizing for incontinence products.



Sizing Tape Measure: Small - XLarge



Sizing Tape Measure: 2XL & Bariatric

- How do I know which size brief will fit me?

Visit: <http://www.prevail.com/x192.php>

First Quality Products is the only manufacturer in the industry with six sizes of adult briefs! Finding your correct size is as easy as 1-2-3.

1. Measure the hips or waist, whatever is larger. Measure hipbone to hipbone and over the stomach.
2. After obtaining the hipbone to hipbone measurement, double the measurement and add two inches.

Example:

Hipbone to hipbone = 22"  
 $22" \times 2 = 44" + 2" = 46"$   
46" measures into a large  
brief

### 3. Sizes of briefs

Youth = 15" - 22" (White)  
Small = 20" - 31" (Green)  
Medium = 32" - 44" (White)  
Regular = 40" - 48" (Lavender)  
Large = 45" - 58" (Blue)  
X-Large = 59" - 64" (Beige)

- How do I put on a First Quality Brief?

AGAIN IT IS AS EASY AS 1-2-3!

### 1. Brief Positioning

Position the brief so that the top tape tabs are lined up with your hipbones.

### 2. Tug Technique

Fold the wings of the front panel back to bring the leg cuffs outward. Pull the panel up snug between the legs until the leg elastic is resting in the crease of the groin area. This will prevent leaking and redness in the inner groin area.

### 3. Tape Tabs

Always fasten the bottom tape tabs first on an upward angle. This will form a secure seal around the buttocks. Fasten the top tabs on a downward angle.