

Specialty Stimulating Electrodes
Electrodes stimulantes de spécialité
Stimulationselektroden für Spezialanwendungen
Elettrodi stimolanti speciali
Electrodos de estimulación para aplicaciones especiales
Stimulerande specialelektroder
Speciale stimulatíe-elektroden
Électrodos de estimulação especiais
Stimuloivat erikoiselektrodit
Specialstimuleringsselektroder
Εξειδικευμένα ηλεκτρόδια διέγερσης
Specializované stimulačné elektrody
Speciális stimuláló elektródok
Специальные стимулирующие электроды
Electrody drażniące do specjalnych zastosowań
Özel Uyarıcı Elektrotlar
Spezielle stimuleringsselektroder
Speciálne stimulačné elektródy

Reusable, Self-adhering for TENS-NMES-FES
Réutilisables adhésives pour TENS/NMES/FES
Wiederverwendbare, selbsthaftende Elektroden für TENS-NMES-FES
Autoadesivi riutilizzabili per TENS-NMES-FES
Autoadhesivo reutilizable para TENS-NMES-FES
Återanvändbar, självhäftande för TENS-NMES-FES
Herbruikbare, zelfklevende elektroden voor TENS-NMES-FES
Auto-adesivos reutilizáveis para TENS/NMES/FES
Toistokäyttöinen, itseliimautuva TENS-NMES-FES:lle
Genanvendelig selvklebende til TENS-NMES-FES
Επαναχρησιμοποίηση, αυτοκόλλητα για TENS-NMES-FES
Samolepící pro více použití pro TENS-NMES-FES
Újrafelhasználható, öntapadó, TENS-NMES-FES kezeléshez
Многоразовые самоклеящиеся для TENS-NMES-FES
Do wielokrotnego użytku, samoprzywierające, do stosowania w zabiegach TENS-NMES-FES
Tekrar kullanılabilir, TENS-NMES-FES için kendiliğinden yapışkan
Selvklebende til flergangsbruk for TENS-NMES-FES
Opakovane použiteľné samolepiace na TENS-NMES-FES

Instructions For Application of Electrodes:

1. Make certain stimulation device is off prior to applying or removing electrodes.
2. Apply electrode only to unbroken skin. Clean electrode application area with soap and water, rinse and dry.
3. With electrodes still on liner, connect lead wire from stimulation device to electrode connector.
4. Remove electrode from liner by grasping the edge of the electrode and peeling it off the liner. Retain liner for storage.
5. Place the electrode on exact skin location (should be determined by a medical professional) by applying the center of the electrode first and smoothing down to electrode edges. Attach lead wire to stimulation device.

Removal of Electrodes:

1. Turn off stimulation device. Disconnect lead wires from stimulation device.
2. Remove electrode from skin by peeling from electrode edge.

Re-Application And Storage of Electrodes:

If adhesive becomes over-saturated, allow electrode to air dry in a refrigerator with adhesive side up until gel regains its tackiness. If the electrode gel appears dry, add a few drops of water to the electrode gel. Let rest to tacky stage and apply to skin. Between uses return the electrode to liner and store in the resealable bag in a cool place out of direct sunlight.

Note: The life of the electrode varies depending on skin conditions, skin preparation, type of stimulation, storage, and climate.

Caution: If skin irritation develops, discontinue use and consult your physician. For single patient use only. These electrodes may be repositioned several times on the same patient.



Caution, consult accompanying documents

Single patient use
Usage limité à un seul patient
Zum Gebrauch an einem einzelnen Patienten
Per uso monopaziente
Para uso en un solo paciente
För enpatientsbruk
Uitsluitend voor gebruik bij één patient
Utilização num único doente
Yhdelle potilaalle

Må kun bruges til én patient
Για χρήση σε έναν ασθενή
Pro jedného pacienta
Egyetlen beteghez történő használatra
Для применения у одного пациента
Do użytku u jednego pacjenta
Tek hastada kullanılmık
Til bruk på én pasient
Na použitie len u jedného pacienta

© 2012 Covidien. Made In USA.

Covidien llc, 15 Hampshire Street, Mansfield, MA 02048 USA.

CEP02 Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.

WB46P0015600