


INSTRUCTION OF USE

Acti-Lance is very easy to handle and its intuitive top button activation minimizes training required for new users



Twist off protective cap and pull it straight out



Place Acti-Lance on the desired puncture site and press the top push button to activate the device.



Gently apply intermittent pressure near the puncture site to obtain the required blood volume.

Recommendations for blood collection test:

- Patients should wash their hands with warm water and soap before the test. Remember that warm fingers increase the blood flow.
- Fingertip test site should be selected slightly off-center to minimize any pain and maximize blood flow.
- The patient's arm should hang down below the elbow to allow blood flow to fingertips.
- After the blood collection, use a tissue and apply firm pressure directly over the finger-stick site.