


Acti-Lance®

Comfort & Safety

SINGLE-USE MEDICAL DEVICES STERILE SAFETY LANCETS

designed for capillary blood sampling
from the fingertip



Product intended use

Acti-Lance® are sterile, single-use devices designed for the sampling of capillary blood from the fingertip of patients. Intended users are healthcare professionals in secondary care and primary care such as doctor's offices. The device is also intended to be used by non-clinical users - patients.

Product Features

Safety

Comfort



Sharps injuries expose healthcare workers to blood-borne pathogens. The Centers for Disease Control and Prevention (CDC) estimates that about 385,000 sharps-related injuries occur annually among healthcare workers in hospitals.

It has been estimated that about half or more of sharps injuries go unreported. Most reported sharps injuries involve nurses, but laboratory workers, physicians, housekeepers, medical students and other healthcare workers are also exposed to injuries.

Sharp injuries can occur before, during or even after use of a sharp and continue to be a serious health and safety hazard.

Sharps injury - facts

Over 80% of needlestick injuries can be prevented with the use of safe needle devices, which, in conjunction with worker education and work practice controls, can reduce injuries by over 90%!*


HTL-STREFA is committed to manufacturing safety products designed to reduce or eliminate the risk of exposure to blood-borne pathogens.

Safety lancets, including **Acti-Lance®**, are an integral component to the sharps injury prevention programs in hospitals, clinics, laboratories, doctor's offices and wherever both patients and professionals need to feel safe.


By using safety lancets, needlestick injuries and infections can be minimized.

**The WHO estimates based on CDC Report*


How to use Acti-Lance®:



Twist off the protective cap and pull it straight out.



Press the **Acti-Lance®** body firmly against the puncture site and press the push-button to activate the device.






Gently apply intermittent pressure near the puncture site to obtain the required blood volume.

Acti-Lance®

Comfort & Safety

Acti-Lance® single-use, sterile safety lancets are push-button activated medical devices used to perform finger pricks in order to collect micro samples of capillary blood.

Available in three versions:

Product version (needle gauge/blade)	Penetration depth
 Acti-Lance® Lite 28G needle	1.5 mm
 Acti-Lance® Universal 23G needle	1.8 mm
 Acti-Lance® Special blade 1.5 mm	2.0 mm

References:

Directions for safety lancets development based on nurses' needs; Poland, June 2016, HTL-STREFA S.A. data on file
<https://www.osha.gov/needlesticks/needlefaq.html>
Mary Foley, MS, RN, Annemarie T. Leyden, EdD, RN American Nurses Association – Independent Study Module
Needlestick Safety and Prevention; http://www.who.int/occupational_health/activities/1anaism.pdf

Poland



HTL-STREFA S.A.
ul. Adamówek 7
95-035 Ozorków
info@htl-strefa.pl
www.htl-strefa.com



Find us on Facebook
www.facebook.com/HTLSTREFA

USA

HTL-STREFA, Inc.
3005 Chastain Meadows Pkwy
Suite 300
Marietta, GA 30066

+1 770 528 0410
info@htl-strefa.com
www.htl-strefa.com



610/ALL/11/2017
Made in Poland

CE mark refers to the product
not to the brochure

CE 0344