

[Print](#) [Close](#)

REPLICARE◇ THIN

Thin Hydrocolloid Dressing

[Overview](#)[Instructions for Use](#)

Description

REPLICARE[®] Thin is a hydrocolloid dressing

Features

REPLICARE[®] Thin is made from a polyurethane film with a thin layer of absorbent colloid.

Benefits

Moist wound environment

REPLICARE[®] Thin maintains a moist wound environment that assists in promoting autolytic debridement while managing low levels of exudate.

Goes on easy, comes off clean

The sterile, no-touch, self-adhering application system makes REPLICARE[®] Thin simple to apply – even while wearing gloves. High cohesive properties keep the dressing from falling apart and leaving residue in the wound, so there's minimal risk of local irritation.

Superior conformability

Like original REPLICARE[®], REPLICARE[®] Thin stretches easily and returns to its original shape more completely than other hydrocolloid dressings – even over joints. This suppleness and flexibility allows the dressings to give without migrating, minimizing granulation tissue injury and damage to periwound skin integrity.

Improves patient comfort

REPLICARE[®] Thin permits increased ease of patient mobility without uncomfortable pulling of the skin, especially over joints. The dressing's thin profile also makes it cosmetically acceptable

Cost-effective choice

REPLICARE[®] Thin promotes more cost-effective use of material and resources and nursing services. The dressing is totally transparent, so there is no need to disturb the dressing to check for saturation or healing progress. The waterproof polyurethane film cover material allows the dressing to be easily cleaned and permits patients to bathe or shower without the need for a dressing change. The dressing also aids in the prevention of bacterial contamination. The no-touch application system results in fewer wasted dressings, and the combination of conformability, adhesion and absorption extends wear time and reduces the risk of wound contamination.

2010 SMITH & NEPHEW
