



### **Stomahesive® PROTECTIVE POWDER** REF25510

This preparation helps form a protective barrier against excoriating discharge on weeping areas.

**DIRECTIONS:** Cleanse the area. Apply a light dusting of powder and remove any excess. Apply Stomahesive® wafer or other product. Powder may also be used to fill spaces between the stoma and skin barrier.

For more information on ConvaTec products, contact your health care professional, or call 1-800-422-8811 (U.S. only). In Canada only, call 1-800-465-6302.

© 2009 ConvaTec Inc.  
Made in UK



ConvaTec Inc.  
Skillman, NJ 08558  
1-800-422-8811