

MEDIHONEY™ HONEYCOLLOID™  
DRESSING WITH ACTIVE *LEPTOSPERMUM*  
HONEY is an absorbent primary wound  
dressing made of Active *Leptospermum*  
Honey (+15) in combination with a hydrocol-  
loidal gelling agent. The hydrocolloidal agent  
is contained within the dressing mass. The  
dressing contains 80%, by weight, Active  
*Leptospermum* Honey. As wound exudate  
is absorbed, the dressing – which has a low  
pH and high osmolarity – forms a gel. This  
assists in maintaining a moist environment  
for optimal wound healing, aids in autolytic  
debridement, and allows for non-traumatic  
removal. Similar to hydrocolloidal pastes, this  
product is not occlusive and does require a  
secondary dressing.