

# Fitting Guide Conveen® Optima & Conveen® Active

Like anything new, using a urisheath and urine bag for the first time or changing to a different type of collecting system can take a bit of getting used to.

Maximize security with the perfect fit – follow these simple steps:



1. Measure the diameter of the penis (if in doubt, choose the smaller diameter)



2. Fit the urisheath leaving a small gap between the tip of the penis and the outlet



3. Holding the tip in place, pull the loop to unroll the urisheath over the penis



4. Gently squeeze the unrolled urisheath around the penis to ensure a secure fit



5. Place the bag on the inner right thigh with the inlet facing inward. Wrap around the leg and fasten



6. Connect to the urisheath by pushing the inlet connector into the outlet on the urisheath

Please note that a urisheath should never be re-used or worn longer than 24 hours.

The urine bag should be emptied regularly when it is approximately 2/3 full.