

[Print](#) [Close](#)

IODOSORB[®] GEL

Cadexomer Gel with Iodine



IODOSORB[®] removes barriers to healing by its triple action antimicrobial, absorbing and desloughing properties

[read more](#)

Overview



Instructions for
Use

Step 1

Open tube and squeeze gel on to a non-adherent gauze to the shape of the ulcer and to a depth of 3 to 4 mm. Position the gel face of the gauze on the wound surface and apply mild pressure to fix in place.

Step 2

Gently take off the existing dressing which should be moistened with saline if it is tending to adhere to the wound surface. Remove any remaining IODOSORB[®] with a gentle stream of saline or sterile water if applicable. Gently blot away any excess fluid before re-applying IODOSORB[®].

Please see Overview Tab for Precautions

2010 SMITH & NEPHEW
