

[Print](#) [Close](#)

ALLEVYN◇ GENTLE BORDER

Gel Adhesive Hydrocellular Foam Dressing with Border



**ALLEVYN◇, The No.1 Foam
Dressing**

The ALLEVYN◇ Range

**Overview**Instructions for
use

Description

ALLEVYN◇ Gentle Border Gel Adhesive Hydrocellular Foam Dressing with Border

Indications

- Chronic and acute exudative wounds
- Full and partial thickness wounds such as pressure ulcers
- Leg ulcers
- Diabetic foot ulcers
- Infected wounds
- Malignant wounds
- Surgical wounds
- First and second degree burns
- Donor sites
- Skin tears
- Fungating ulcers

ALLEVYN◇ Gentle Border is suitable for use on fragile skin.

Precautions

- Do not use ALLEVYN◇ Gentle Border dressing with oxidizing agents such as hypochlorite solutions (e.g. EUSOL) or hydrogen peroxide, as these can break down the absorbent polyurethane component of the dressing.
- If reddening or sensitization occur discontinue use.

Features

2010 SMITH & NEPHEW

- Advanced triple-layered construction
- Absorbent hydrocellular pad
- Perforated silicone gel adhesive wound contact layer
- Adhesive border for gentle but effective fixation
- Highly permeable waterproof outer film