

560 Urethral Catheter Tray Instructions and Ingredients

Aplicare Three Benzalkonium Chloride Swabsticks

Product Information
Use Skin cleanser
Warnings For external use only Do not use in the eyes Ask a doctor before use if injuries are: <ul style="list-style-type: none">• Deep or puncture wounds• Serious burns Stop use and ask a doctor if: <ul style="list-style-type: none">• Infection occurs• Redness irritation, swelling, or pain persists or increases Keep out of reach of children. In case of accidental ingestion, seek professional assistance or consult a poison control center immediately
Directions Apply locally as needed
Other information Latex-free
Ingredients Benzalkonium chloride, citric acid anhydrous, disodium hydrogen phosphate anhydrous, water

Aplicare Five Gram Lubricating Jelly

Sterile unless opened or damaged
For professional and hospital use only
Use: For lubrication during catheter insertion
Water soluble