


Hydrocolloid Wound Dressing

ZeniPlast is a sterile, hydrocolloid wound dressing used to help maintain an optimal moist wound environment conducive to wound healing. ZeniPlast interacts with wound exudate to form a soft gel which helps to protect and prevent damage to newly formed tissue during dressing removal.



Product Features

- Helps to create a controlled wound environment
- Low-profile edge helps to prevent lifting
- Extended wear dressing (up to 7days)
- Protective, occlusive barrier
- Low friction cover
- Flexible and conformable for hard to dress areas
- Waterproof
- Gel formation helps to protect fragile tissue
- Easy application

Indications for Use

- Moist to Low Draining wounds
- Dermal Ulcers
- Superficial Wounds
- Post-operative wounds
- Surgical Wounds
- Donor Sites
- Lacerations and Abrasions
- 1st and 2nd Degree Burns

Contraindications

ZeniPlast is contraindicated for use on 3rd Degree burns

Item Number	Description	Size	Packaging	HCPCS Code
50022	ZeniPlast	2"x 2"	10/Box	A6237
50044	ZeniPlast	4"x 4"	10/Box	A6237
50066	ZeniPlast	6"x 6"	10/Box	A6238

*Providers are responsible for determining appropriate billing codes for Medicare Program Claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.