

ZeniFoam-Ag

Foam Wound Dressings with Silver

ZeniMedical



Foam Dressings with Silver

ZeniFoam-Ag dressings are highly absorbent foam dressings with 0.35-0.4mg/cm³ ionic silver to provide an optimal moist wound healing environment combined with effective antibacterial activity and exudate management.

Product Features



- Highly Absorbent
- Can be Cut to Shape
- Primary or Secondary Dressing
- Maintains Moist Environment
- Broad spectrum Antimicrobial against Gram negative and Gram positive bacteria
- Non-Adherent to Wound Bed
- Protects Against Maceration
- Easy Application
- Bacterial reduction up to 7 days

Indications for Use

- Partial Thickness Wounds
- Full Thickness Wounds
- Leg Ulcers
- Pressure Ulcers
- Diabetic Foot Ulcers
- Second Degree Burns
- Donor Sites
- Skin Abrasions
- Post-Op Wounds

Contraindications

ZeniFoam-Ag foam dressings with silver are not intended for use on third-degree burns and should not be used on patients with known sensitivity to silver.

Item Number	Description	Size	Packaging	HCPCS Code
10022-Ag	ZeniFoam-Ag Silver Foam No Adhesive/No Border	2" x 2"	10/Box	A6209
10045-Ag	ZeniFoam-Ag Silver Foam No Adhesive/No Border	4" x 5"	10/Box	A6210
10066-Ag	ZeniFoam-Ag Silver Foam No Adhesive/No Border	6" x 6"	10/Box	A6210
10088-Ag	ZeniFoam-Ag Silver Foam No Adhesive/No Border	8" x 8"	10/Box	A6211

* Providers are responsible for determining appropriate billing codes for Medicare Program Claims. This information is included for convenience and is not intended, no should be considered billing advice.