


Silicone Wound Contact Dressing

ZeniContact Silicone Wound Contact Dressing is an open mesh, soft silicone wound contact layer that allows exudate to pass through to the secondary dressing. It is designed for the management of dry to heavily exuding wounds including venous and arterial ulcers, first and second degree burns, radiation injuries, abrasions, skin tears, skin grafts, blistering, surgical incisions, diabetic ulcers and lacerations.

Product Features



- Conforms Easily
- Flexible and Soft
- 100% Open Mesh
- Allows Exudate Pass-through
- Gentle, Non-adherent Contact Layer
- Can be Cut to Shape
- Primary Dressing
- Prevents Secondary Dressings from Adhering to Wound
- Easy Peel of Release Liner
- Removable Protective Backing Prevents Dressing from Sticking to Gloves During Application
- Can be Used with NPWT Dressings or as Topical Ointment Carrier

Indications for Use

- Traumatic Wounds
- Abrasions
- Blisters
- 1st & 2nd Degree Burns
- Surgical Wounds
- Lacerations
- Pressure Ulcers.
- Skin Tears
- Skin Grafts
- Venous Ulcers
- Radiation Wounds
- Granulation Tissue
- Diabetic Ulcers
- Arterial Ulcers

Contraindications

ZeniContact is not intended for use on 3rd Degree Burns

Item Number	Description	Size	Packaging	HCPCS Code
70047	ZeniContact: 4 x 7	4" x 7"	10/Box	A6207

* Providers are responsible for determining appropriate billing codes for Medicare Program Claims. This information is included for convenience and is not intended, no should be considered billing advice.