Pro-Stat® Sugar Free AWC is a ready-todrink medical food providing 17 g of enzymehydrolyzed complete protein and 100 calories per 1 fl oz.

### Main Features

- Clinically shown in one clinical study to promote pressure ulcer and wound healing\*
- Contains 17 g of protein and 100 calories in 1 oz with added Arginine, Citrulline, Cystine, Vitamin C, and Zinc
- Hydrolyzed collagen protein: high in conditionally essential amino acids
- Sugar free

# **Typical Pro-Stat Amino Acid Profile**

(g/100 g of protein)

Amino acid profile in collagen protein can naturally vary.

## **ESSENTIAL AMINO ACIDS**

| Histidine                 | 0.53  |
|---------------------------|-------|
| Isoleucine                | 1.33  |
| Leucine                   | 2.71  |
| Lysine                    | 3.36  |
| Methionine                | 0.80  |
| Phenylalanine             | 1.86  |
| Threonine                 | 1.86  |
| Tryptophan                | 0.46  |
| Valine                    | 2.14  |
| NON-ESSENTIAL AMINO ACIDS |       |
| Alanine                   | 7.89  |
| Aspartic acid             | 4.63  |
| Serine                    | 2.85  |
| CONDITIONALLY-ESSENTIAL   |       |
| AMINO ACIDS               |       |
| Arginine                  | 20.05 |
| Cystine                   | 0.23  |
| Glutamic acid             | 7.32  |
| Glycine                   | 19.58 |
| Proline                   | 11.14 |
| Tyrosine                  | 0.37  |
| MODIFIED AMINO ACIDS      |       |
| Hydroxylysine             | 0.65  |
| Hydroxyproline            | 10.22 |

## Storage:

- Store at room temperature
- · Refrigeration causes gelling
- Discard 3 months after opening
- Record date opened on bottom of container
- Do Not Use If Seal Is Broken

For best by date and lot number, see bottom of bottle.




# NUTRICIA

# Pro-Stat<sup>®</sup> Sugar Free AWC

# **Advanced Wound Care**

Complete Liquid Protein with Arginine, Citrulline, Cystine, Vitamin C and Zinc





Net content 30 fl oz (887 mL)

Lee S.K., et al. Adv Skin Wound Care 2006; 19:92-96. Pro-Stat is a medical food. Contains 0% Fruit Juice

# Contains 0% Fruit Juice

Pro-Stat® Sugar Free AWC is a medical food designed for the dietary management of wounds, and other conditions requiring increased protein in low volume. Use under medical supervision.

| Nutrients           | per 30 mL serving |
|---------------------|-------------------|
| Calories            | 100 kcal          |
| Calories per mL     | 3.3 kcal          |
| Protein             | 17 g              |
| Includes L-Arginine | 3.3 g             |
| Fat                 | 0 g               |
| Carbohydrate        | 7 g               |
| Vitamin C           | 175 mg            |
| Zinc                | 10 mg             |
| Sodium              | 50 mg             |
| Potassium           | 20 mg             |
| Phosphorus          | 50 mg             |

Ingredients: Hydrolyzed Collagen Protein, Water, Glycerin, Malic Acid, L-Arginine. Contains 2% or less: Natural Flavors, L-Citrulline, L-Tryptophan, L-Cysteine HCL, Phosphoric Acid, Ascorbic Acid, Potassium Sorbate (to help protect flavor), Sodium Benzoate (to help protect flavor), Sucralose, Zinc Sulfate.

Distributed by: Nutricia North America, Rockville, MD 20850 1.800.221.0308 • www.pro-stat.com

# Lactose Free, Soy Free & Gluten Free

# **Viscosity: Honey like**

Indications: For stage 3 & 4 pressure ulcers, hard to heal wounds (diabetic, venous, surgical, burns), and other conditions requiring increased protein needs in low volume such as multiple pressure ulcers, unintentional weight loss, protein-energy malnutrition, low serum proteins, & muscle loss (sarcopenia).

Administration: Use under medical supervision. Healthcare professionals should determine intake based on condition, estimated need, and dietary intake. Suggested intake: 1-3 servings per day. Administer orally or through feeding tube. Not suitable for use as a sole source of nutrition. Not for parenteral use. Suitable for adults.

Oral Feeding: Can be taken straight or mixed with any hot or cold beverages or food.

# **Tube Feeding:**

- 1. Flush feeding tube with 30-60 mL of water.
- 2. Pour 30 mL of Pro-Stat in a 4-6 fl oz container.
- 3. Add 30-60 mL water and mix well with disposable spoon or tongue blade.
- 4. Administer Pro-Stat via syringe.
- 5. Flush with 30-60 mL of water.

Manufactured under U.S. Patents Nos. 7,318,934, 7.682,624, & 7.682,625.

Item # 40130 Reimbursement Code 26974-0410-31 HCPCS Code B4155

