




# ZeniFoam and ZeniFoam Gentle

Polyurethane Foam Wound Dressings

## Polyurethane Foam Wound Dressings

ZeniFoam and ZeniFoam Gentle are sterile polyurethane foam dressings available with or without Silicone-Faced adhesives and borders. ZeniFoam is highly absorbent and helps to protect sensitive and fragile skin.



## Product Features

- Highly Absorbent
- Non-adherent to wound bed
- Lint free, no particulates left in wound
- Available in multiple sizes and styles
- Helps protect against maceration
- Gentle on sensitive and fragile skin
- Highly conformable and comfortable
- Can be repositioned upon application
- Water resistant backings
- Easy application
- Can be used as Primary or Secondary Dressing

## Indications for Use

- Pressure Ulcers
- Partial Wounds
- Full-Thickness Wounds
- Leg Ulcers
- Diabetic Ulcers
- Surgical Wounds
- Donor Sites
- Lacerations and Abrasions
- 1<sup>st</sup> and 2<sup>nd</sup> Degree Burns
- Moderate to Highly draining wounds

## Contraindications

ZeniFoam is contraindicated for use on 3<sup>rd</sup> Degree Burns or as an internal dressing/sponge

Item Number	Description	Size	Packaging	HCPCS Code
10022	ZeniFoam (No Border, No Adhesive)	2"x 2"	10/Box	A6209
10045	ZeniFoam (No Border, No Adhesive)	4"x 5"	10/Box	A6210
10066	ZeniFoam (No Border, No Adhesive)	6"x 6"	10/Box	A6210
20022	ZeniFoam Gentle (Silicone-Faced)	2"x 2"	10/Box	A6209
20044	ZeniFoam Gentle (Silicone-Faced)	4"x 4"	10/Box	A6209
20066	ZeniFoam Gentle (Silicone-Faced)	6"x 6"	10/Box	A6210
30022	ZeniFoam, Gentle Border (Silicone)	2"x 2"	10/Box	A6412
30044	ZeniFoam, Gentle Border (Silicone)	4"x 4"	10/Box	A6212
30066	ZeniFoam, Gentle Border (Silicone)	6"x 6"	10/Box	A6213

\*Providers are responsible for determining appropriate billing codes for Medicare Program Claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.