

ZeniFoam Dressings

Polyurethane Foam Dressings

Available in 3 different styles: Foam without adhesive border; Foam with silicone non-adherent contact layer without border; or Foam with silicone non-adherent contact layer with a breathable, conformable adhesive coated fabric tape.

STERILE

Single Patient Use Only, Disposable

Not made with Natural Rubber Latex



Ingredients:

ZeniFoam:

Polyurethane Foam dressing contains CMC (carboxymethylcellulose) with waterproof backing made of oxygen permeable polyurethane film.

ZeniFoam Gentle:

Polyurethane Foam dressing contains CMC (carboxymethylcellulose) with non-adherent Silicone contact layer. Waterproof backing made with oxygen permeable polyurethane film.

ZeniFoam Gentle Border:

Polyurethane Foam dressing contains CMC (carboxymethylcellulose) with nonadherent Silicone contact layer. Waterproof backing made with oxygen permeable polyurethane film. Tape border consists of gentle adhesive on woven fabric.

Manufactured for:

ZeniMedical

1410 East Weisgarber Road

Knoxville, Tennessee 37909

Part #: 100200300

Rev. 10/17