
ZeniFiber

Calcium Alginate Wound Dressing



Calcium Alginate

ZeniFiber is a sterile, nonwoven pad or rope made from calcium alginate fibers. Upon contact with wound exudate, ZeniFiber becomes a fibrous gel that helps to maintain a moist wound environment that is conducive to wound healing.



Product Features

- Highly absorbent
- Excellent fluid handling capacity
- Easy dressing changes
- Helps to maintain moist wound environment
- Can be used on deep or shallow wounds
- Used as a primary dressing
- For use on Moderate to Heavily draining wounds

Indications for Use

- Pressure Ulcers
- Surgical Wounds

Donor Sites

- Partial Thickness Wounds•
 - Full-Thickness Wounds Lacerations and Abrasions
- Leg Ulcers

- 1st and 2nd Degree Burns
- Diabetic Ulcers

Contraindications

ZeniFiber is contraindicated for use on 3rd Degree Burns, as an internal would Dressing or sponge, or for Patients with known sensitivity to alginates

<u>Item Number</u>	Description	<u>Size</u>	<u>Packaging</u>	HCPCS Code*
60022	ZeniFiber	2"x 2"	10/Box	A6196
60045	ZeniFiber	4"x 5"	10/Box	A6197
60066	ZenfFiber	6" Rope	10/Box	A6199

^{*}Providers are responsible for determining appropriate billing codes for Medicare Program Claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.